



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Mobile Shadowless Operation Lamp Single Dome Stand
Model With Castors & With Battery Backup (4 Bulbs)

Product Code :

SAMCOMEDICAL200001A074

**Description :**

Mobile Shadowless Operation Lamp Single Dome Stand Model With Castors & With Battery Backup (4 Bulbs)
Illumination Intensity : 40000 Lux Color Temperature : 3500-5000 K Spot Diameter : 100-200mm Height
Adjustment : 400mm Lamphead Diameter : 500mm Power Supply : 220V / 50Hz Power Input : 120 VA With
battery backup



Samco Medical

Website: www.samcomedical.com,