



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Mobile Shadowless Operation Lamp Single Dome Stand
Model With Castors (5 Bulbs)

Product Code :

SAMCOMEDICAL200001A075

**Description :**

Mobile Shadowless Operation Lamp Single Dome Stand Model With Castors (5 Bulbs) Illumination Intensity : 50000 Lux Color Temperature : 3500-5000 K Spot Diameter : 100-200 mm Height Adjustment : 400 mm Lamphead Diameter : 500 mm Power Supply : 220V / 50Hz Power Input : 150 VA



Samco Medical

Website: www.samcomedical.com,