



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
LED Surgical Light with Castors Floor Model (For Minor Surgical Procedures)

Product Code :
SAMCOMEDICAL200001A079



Description :

LED Surgical Light with Castors Floor Model (For Minor Surgical Procedures) Illumination Intensity : 40000 Lux
Color Temperature : 6000-6500 K Color rendering index : 75 (Ra) Spot diameter : 150mm Light Field Diameter :
D 30 cm = 18 cm Temperature Rise of Field : 0 ($1\frac{1}{2}^{\circ}\text{C}$) Power Supply : 220 V / 50 Hz. AC Extention : 75 cm
goose neck Specification of Bulb : LED 6 bulbs (high power) Brightness Control : Available



Samco Medical

Website: www.samcomedical.com,