



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

LED Examination Light With Stand Floor Model & With
Microcomputer Controlled 8 Step Brightness Adjustment

Product Code :

SAMCOMEDICAL200001A083

**Description :**

LED Examination Light With Stand Floor Model & With Microcomputer Controlled 8 Step Brightness Adjustment
Standard Features : Microcomputer control 8-steps brightness adjustment Technical specifications Color
Temperature : 3000-5500 K Illumination : ? 20000 Lux Temperature Rise of Field : ?8($i\frac{1}{2}$ C) Power Supply : 220
V / 50 Hz. AC Power (W) : ?12 Type of bulb : Super bright LED Luminotron 3V/0.5Wx19 Spot Diameter : 80-120
mm Range : 300-500mm from lamp adapter Brightness Control : 8 step digital control



Samco Medical

Website: www.samcomedical.com,