



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

LED Examination Light with Stand Floor Model &
Microcomputer Controlled 8 Step Brightness Adjustment
With Battery Backup

Product Code :

SAMCOMEDICAL200001A084

**Description :**

LED Examination Light with Stand Floor Model & Microcomputer Controlled 8 Step Brightness Adjustment With Battery Backup Standard Features: Microcomputer control 8-steps brightness adjustment Technical specifications Color Temperature : 3000-5500 K Illumination : ? 20000 Lux Temperature Rise of Field : ?8(1/2C) Power Supply : 220 V / 50 Hz. AC Power (W) : ?12 Type of bulb : Super bright LED Luminotron 3V/0.5Wx19 Spot Diameter : 80-120 mm Range : 300-500mm from lamp adapter Brightness Control : 8 step digital control



Samco Medical

Website: www.samcomedical.com,