



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
LED Examination Light Floor Model With Castors

Product Code :
SAMCOMEDICAL200001A086



Description :

LED Examination Light Floor Model With Castors Color Temperature : 3000-5500 K Illumination : ? 3000 Lux
Temperature Rise of Field : ? 8 ($1\frac{1}{2}$ C) Power Supply : 220 v / 50 Hz. AC Power : ? 1.2 (W) Type of bulb : Super
bright LED Luminotron 3V/0.2W x 8 Spot Diameter : 80-120mm Range : 300-500mm from lamp adapter
Brightness Control : Switch ON/OFF



Samco Medical

Website: www.samcomedical.com,