



Email : sales@samcomedical.com  
Phone: 1-300-400-8211

**Product Name :**  
LED Examination Light Table Top Model

**Product Code :**  
SAMCOMEDICAL200001A087



**Description :**

LED Examination Light Table Top Model Technical specifications Color Temperature : 3000-5500 K Illumination : ? 12000 Lux Temperature Rise of Field : ? 8 (C) Power Supply : A.C 220V 50Hz Power : ? 3.5 (W) Type of bulb : Super bright LED Luminotron 3V/0.5Wx8 Spot Diameter : 80-120mm Range : 300-500mm from lamp adapter Brightness Control : Digital



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),