



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Suction Unit Eurovac AC/DC

Product Code :
SAMCOMEDICAL150001A003



Description :

Suction Unit Eurovac AC/DC

Technical Specification :

Suction Unit Eurovac AC/DC Housing: ABS Plastic Moulded Capacity: -600 mm Hg 10 , 55 Ltrs/ min Pump Type: Diaphragm Pump Jars: 1 x 2 Ltrs Polycarbonate Jar Filter: Bacterial Filter Autoclavable/Reusable Tubing: 8 mm ID x 2 mtr. (Silicon) Vacuum Gauge: 6.25 cm dia 0-760 mm Hg calibration. Power: 220 V AC, 50 / 60 Hz, 70Watt, 12 V DC (Internal Rechargeable Battery 12V 4.5Amp) Battery: Lead Acid (Ni-cd Optional) Noise Level: 55 dB A 3 Dimension & Weight: 43 x 22 x 37 cms, 11 Kg. Overflow Safety: Mechanical / Electronic Acoustic Alarm Accessories: Car Charger & Power Adaptor



Samco Medical

Website: www.samcomedical.com,