



Email : sales@samcomedical.com  
Phone: 1-300-400-8211

**Product Name :**  
Suction Unit Eurovac AC

**Product Code :**  
SAMCOMEDICAL150001A007



**Description :**

Suction Unit Eurovac AC

**Technical Specification :**

Suction Unit Eurovac AC Housing: ABS Plastic Moulded Capacity: -575 mm Hg 10 , 45 Ltrs/ min Pump Type: Diaphragm Pump Jars: 1 x 2 Ltrs Polycarbonate Jar Filter: Bacterial Filter Autoclavable/Reusable Tubing: 8 mm ID x 2 mtr. (Silicon) Vacuum Gauge: 6.25 cm dia 0-760 mm Hg calibration. Power: 220 V AC, 50 / 60 Hz, 90Watt (110V on request) Noise Level: 55 dB A 3 Dimension & Weight: 43 x 22 x 37 cms, 8 Kg. Overflow Safety: Mechanical / Electronic Acoustic Alarm



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),