



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Samcomedical Manual Suction Unit Foot with Pedal
Operated

Product Code :
SAMCOMEDICAL150005A003



Description :

Samcomedical Manual Suction Unit Foot with Pedal Operated

Technical Specification :

Manual Suction Unit (Foot / Pedal Operated) Capacity: Foot Operated Vacuum pump giving vacuum more than 550mm (+ 50mm). Jars: Two Autoclavable PC jars with cap of 1 litre capacity. A suction cannula is also provided with the unit. The unit is fitted on a powder coated base



Samco Medical

Website: www.samcomedical.com,