



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Suction Unit - Dual Power (Electric cum Manual)

**Product Code :**  
SAMCOMEDICAL150002A001

**Description :**

Suction Unit - Dual Power (Electric cum Manual) Housing: MS Powder Coated, SS Top on two wheels Capacity (Electrical): -600mmHg  $\dot{\text{I}}\dot{\text{C}}\frac{1}{2}$  10 , 45 Ltrs/ min Capacity (Pedal): -600mmHg  $\dot{\text{I}}\dot{\text{C}}\frac{1}{2}$  10 , 200ml/ Stroke Pump Type: Electric - Diaphragm / Pedal - Piston Jars: 2 x 2 Ltrs Polycarbonate Jar Filter: Bacterial Filter Autoclavable/Reusable Tubing: 10mmID x 2 mtr. (PVC) Vacuum Gauge: 6.25 cm dia 0-760mmHg calibration. Power: 220 V AC, 50 / 60 Hz, 90Watt (110V on request) Noise Level: 55 dB A  $\dot{\text{I}}\dot{\text{C}}\frac{1}{2}$  3 Dimension & Weight: 42 x 39 x 64 cms, 15.5 Kg Overflow Safety: Dual Mechanical



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),