



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Suction Unit HI-VAC Plus (90 Ltrs./Min)

Product Code :
SAMCOMEDICAL150004A001

Description :

Suction Unit HI-VAC Plus (90 Ltrs./Min) Housing: Plastic moulded on two wheels Capacity: -700 mm Hg $\pm$ 10, 90 ltrs/ min Pump Type: Double Piston Pump Jars: 2 x 2.5 Ltrs wide mouthed (Polycarbonate) with changeover lever. Provided with handle for hygiene while emptying/cleaning. Filter: Bacterial Filter Autoclavable/Reusable Tubing: 10 mm ID x 2 mtr. (PVC) Vacuum Gauge: 6.25 cm dia 0-760 mm Hg calibration. Power: 220 V AC, 50 / 60 Hz, 580Watt (110V on request) Noise Level: 60 dB A $\pm$ 3 Dimension & Weight: 46 x 34 x 86 cms, 23 Kg. Overflow Safety: Mechanical / Electronic Acoustic Alarm Accessories: Foot Switch



Samco Medical

Website: www.samcomedical.com,