



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Intensive Care Hospital Bed HI-LOW

Product Code :
SAMCOMEDICAL1001A001

Description :

Intensive Care Hospital Bed HI-LOW

Technical Specification :

Size : 198 L x 80 W cm. - Height adjustable from 60 cms. to 80 cms. Four Section steel sheet top uniformly perforated. Three separate screw levers for manoeuring backrest, kneerest and HI-LOW positions. Head and leg ends made of S.S. tubular pipes assembled to laminated panels. For rubber buffers provided to prevent damage to hospital walls. Mounted on four castors (two with brakes). Pre treated and epoxy powder coated.



Samco Medical

Website: www.samcomedical.com,