



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Suction Unit MAGNAVAC (70 Ltr/Min)

Product Code :
SAMCOMEDICAL150004A003



Description :

Suction Unit MAGNAVAC (70 Ltr/Min)

Technical Specification :

Suction Unit MAGNAVAC (70 Ltr/Min) High Vacuum Suction Unit for LIPO SURGERY Housing: Stainless Steel on four castors Capacity: -700 mm Hg 10 , 70 Ltrs/ min Pump Type: Double Rotary Vane Pump Jars: 2 x 2 Ltrs Polycarbonate Jar Filter: Trap Filter Tubing: 10 mm ID x 2 mtr. (PVC) Vacuum Gauge: 10 cm dia 0-760 mm Hg calibration. Power: 220 V AC, 50 / 60 Hz, 380Watt (110V on request) Noise Level: 50 dB A 3 Dimension & Weight: 49 x 36 x 70 cms, 36 Kg. Overflow Safety: Mechanical



Samco Medical

Website: www.samcomedical.com,