



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Suction Unit Mobile

Product Code :
SAMCOMEDICAL150004A005



Description :

Suction Unit Mobile

Technical Specification :

Suction Unit - Mobile Housing: MS Powder Coated, SS Top on two wheels Capacity: -700 mm Hg 10 , 25 Ltrs/ min Pump Type: Single Rotary Vane Pump Jars: 2 x 2 Ltrs Polycarbonate Jar Filter: Trap Filter Tubing: 8 mm ID x 2 mtr. (PVC) Vacuum Gauge: 6.25 cm dia 0-760 mm Hg calibration. Power: 220 V AC, 50 / 60 Hz, 180Watt (110V on request) Noise Level: 55 dB A 3 Dimension & Weight: 43 x 30 x 62 cms, 24 Kg Overflow Safety: Mechanical



Samco Medical

Website: www.samcomedical.com,