



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Enema/Douche Tube & Fitting Set

Product Code :
SAMCOMEDICAL140001A017

Description :

Enema/Douche Tube & Fitting Set

Technical Specification :

Enema/Douche Tube & Fitting Set Hose adapter. Superior white rubber hose (tubing). Pinch clamp.
Douche/vaginal tip. Enema/Rectal tip. Hanging Hook(A multi purpose system with attachments that serves as water bottle, douche or enema syringe)



Samco Medical

Website: www.samcomedical.com,