



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

S.M.A.R.T. Baby Printed Rubber Sheeting

Product Code :

SAMCOMEDICAL140001A029

Description :

S.M.A.R.T. Baby Printed Rubber Sheeting

Technical Specification :

S.M.A.R.T. Baby Printed Rubber Sheeting Made from premium quality technically specified imported rubber Silicone coated for extra smoothness to the baby's tender skin. Autoclavable with steam, hot air, ETO, gamma rays. Washable with luke warm water and mild detergents. Freedom from wrinkles and fold marks. Highly durable



Samco Medical

Website: www.samcomedical.com,