



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Infant Radiant Warmer with Fixed Baby Cradle and Drawers (Standard) with Microprocessor based temperature controller with 3 modes Skin / Air / Manual mode

Product Code :

SAMCOMEDICAL410001A015

**Description :**

Infant Radiant Warmer with Fixed Baby Cradle and Drawers (Standard) with Microprocessor based temperature controller with 3 modes Skin / Air / Manual mode These are mobile Warmer mounted on wheels with fixed baby cradle, lower shelf and 3 drawers storage Cabinet. Mechanical Specification Heater Box Consists of Ceramic Infra Red Heater 650 Watts Metallic overhead heater box. The construction of the stainless steel reflector and Guard rail is such that the complete cleaning of heater & reflector is possible without any electrical shock hazards Halogen Examination Lamp (12 Volts, 50 Watts) To provide accurate assessment of Infant colour and sufficient illumination during procedures. Baby Tray (Stainless Steel) Clear acrylic collapsible side panels . Stainless steel baby tray with transparent slot to incorporate the facility of under surface phototherapy. Acrylic tray and cushion mattress with head up/down facility. Acrylic X-ray tray with X-ray guide for positioning the cassette. The design aspects also allow ease of dismantling baby care area during cleaning or fumigation Essential Attachments Lower Shelf and 3 Drawer Storage Cabinet. S.S. IV Stand



Samco Medical

Website: www.samcomedical.com,