



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Infant Radiant Warmer with Fixed Baby Cradle (Premium)
with Microprocessor based temperature controller with 2
modes Skin / Manual mode

Product Code :

SAMCOMEDICAL410001A018

**Description :**

Infant Radiant Warmer with Fixed Baby Cradle (Premium) with Microprocessor based temperature controller with 2 modes Skin / Manual mode This is mobile Warmer mounted on wheels with fixed baby cradle. Mechanical Specification Heater Box Consists of Ceramic Infra Red Heater 650 Watts The source box can be swiveled to take X-ray. Due to temp. insulated cover on the heater box, its exterior surface remains cool even after long run. The construction of the stainless steel reflector and guard rail is such that the complete cleaning of heater & reflector is possible without any electrical shock hazards Halogen Examination Lamp (12 Volts, 50 Watts) To provide accurate assessment of Infant colour and sufficient illumination during procedures. Baby Tray (Stainless Steel) Clear acrylic collapsible side panels . Stainless steel baby tray with transparent slot to incorporate the facility of under surface phototherapy. Acrylic tray and cushion mattress with head up/down facility. Acrylic X-ray tray with X-ray guide for positioning the cassette. The design aspects also allow ease of dismantling baby care area during cleaning or fumigation Essential Attachments 2 Lower Shelves. Side Mounted Monitor tray, Pole for Syringe Pump / Infusion Pump & S.S. IV Stand



Samco Medical

Website: www.samcomedical.com,