



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**

Radiant Warmer Stand / Radiant Overhead Lamp  
(Premium) with Microprocessor based temperature controller  
with 2 modes Skin / Manual mode

**Product Code :**

SAMCOMEDICAL410001A022

**Description :**

Radiant Warmer Stand / Radiant Overhead Lamp (Premium) with Microprocessor based temperature controller with 2 modes Skin / Manual mode This is mobile Warmer Stand mounted on wheels. Mechanical Specification Heater Box Consists of Ceramic Infra Red Heater 650 Watts The source box can be swivelled to take X-ray. Due to temp. insulated cover on the heater box, its exterior surface remains cool even after long run. The construction of the stainless steel reflector and Guard rail is such that the complete cleaning of heater & reflector is possible without any electrical shock hazards Halogen Examination Lamp (12 Volts, 50Watts) To provide accurate assessment of Infant colour and sufficient illumination during procedures Essential Attachments Side Mounted Monitor tray, Pole for Syringe Pump / Infusion Pump. S.S. IV Stand



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),