



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Infant Radiant Warmer with Detachable Baby Cradle and Drawers (Standard) with Microprocessor based temperature controller with 2 modes Skin / Manual model

Product Code :

SAMCOMEDICAL410001A024

**Description :**

Infant Radiant Warmer with Detachable Baby Cradle and Drawers (Standard) with Microprocessor based temperature controller with 2 modes Skin / Manual model This is mobile Warmer Stand mounted on wheels. Technical Specification Heater Box Consists of Ceramic Infra Red Heater 650Watts. Metallic overhead heater box. The construction of the stainless steel reflector and guard rail is such that the complete cleaning of heater & reflector is possible without any electrical shock hazards Halogen Examination Lamp (12 Volts, 50 Watts) To provide accurate assessment of Infant colour and sufficient illumination during procedures. Infant Care Trolley Standard Clear acrylic collapsible side panels to assess an infant from all sides Stainless steel baby tray with transparent slot to incorporate the facility of under surface phototherapy Acrylic tray and cushion mattress with head up/down facility Acrylic X-ray tray with X-ray guide for positioning the cassette Two Lower shelves Mounted on wheel 2" The design aspects also allow easy of dismantling baby care area during cleaning or fumigation Essential Attachments S.S. IV Stand



Samco Medical

Website: www.samcomedical.com,