



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**

Radiant Warmer Stand / Radiant Overhead Lamp  
(Standard) with Microprocessor based temperature controller  
with 3 modes Skin / Air / Manual model

**Product Code :**

SAMCOMEDICAL410001A025

**Description :**

Radiant Warmer Stand / Radiant Overhead Lamp (Standard) with Microprocessor based temperature controller with 3 modes Skin / Air / Manual model This is mobile Warmer Stand mounted on wheels. Mechanical Specification Heater Box Consists of Ceramic Infra Red Heater 650Watts. Metallic overhead heater box. The construction of the stainless steel reflector and guard rail is such that the complete cleaning of heater & reflector is possible without any electrical shock hazards Halogen Examination Lamp (12 Volts, 50 Watts) To provide accurate assessment of Infant colour and sufficient illumination during procedures Essential Attachments Lower Shelf. S.S. IV Stand



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),