



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**

Infant Radiant Warmer with Fixed Baby Cradle and Drawers with Microprocessor based temperature controller with 3 modes Skin / Air / Manual model

**Product Code :**

SAMCOMEDICAL410001A027

**Description :**

Infant Radiant Warmer with Fixed Baby Cradle and Drawers with Microprocessor based temperature controller with 3 modes Skin / Air / Manual model These are mobile Radiant Warmer fixed height mounted on wheels with fixed baby cradle, lower shelf and 3 drawers storage Cabinet. Technical Specification Heater Box Heating Element : Ceramic Heater 650 W with parabolic reflector, protected by metallic grid The source box can be swivelled to take X-ray Halogen Examination Lamp (12 Volts, 50 Watts) Provides accurate assessment of Infant colour and sufficient illumination during procedures Transparent acrylic, drop down and lockable side panels Baby Tray with mattress, along with head up/down facility and X-ray cassette holder Mattress density approx 21-25Kg/m<sup>3</sup> and with removable, washable, waterproof cover Acrylic X-ray tray with X-ray guide for positioning the cassette The design aspects also allow ease of dismantling baby care area during cleaning or fumigation Essential Attachments Lower Shelf and plastic moulded 3 Storage Drawer Side tray for Monitor Pole for Syringe Pump / Infusion Pump S.S. IV Rod 4 swivel castors (100mm) - 2 castors with brake



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),