



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**

Infant Radiant Warmer with Fixed Baby Cradle with  
Microprocessor based temperature controller with 3 modes  
Skin / Air / Manual

**Product Code :**

SAMCOMEDICAL410001A029

**Description :**

Infant Radiant Warmer with Fixed Baby Cradle with Microprocessor based temperature controller with 3 modes Skin / Air / Manual This is mobile Warmer mounted on wheels with fixed baby cradle. Technical Specification Heater Box Heating Element : Ceramic Heater 650 W with parabolic reflector, protected by metallic grid The source box can be swivelled to take X-ray Halogen Examination Lamp (12 Volts, 50 Watts) Provides accurate assessment of Infant colour and sufficient illumination during procedures Transparent acrylic, drop down and lockable side panels Baby Tray with mattress, along with head up/down facility and X-ray cassette holder Mattress density approx 21-25Kg/m<sup>3</sup> and with removable, washable, waterproof cover Acrylic X-ray tray with X-ray guide for positioning the cassette The design aspects also allow ease of dismantling baby care area during cleaning or fumigation Essential Attachments 2 Lower Shelves Side tray for Monitor Pole for Syringe Pump / Infusion Pump S.S. IV Rod 4 swivel castors (100mm) - 2 castors with brake



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),