



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Infant Radiant Warmer with Detachable Baby Cradle with Microprocessor based temperature controller with 3 modes Skin/Air/Manual

Product Code :

SAMCOMEDICAL410001A031

**Description :**

Infant Radiant Warmer with Detachable Baby Cradle with Microprocessor based temperature controller with 3 modes Skin/Air/Manual This is mobileWarmer Stand mounted on wheels with detachable infant care trolley
Technical Specification Heater Box Heating Element : Ceramic Heater 650Wwith parabolic reflector, protected by metallic grid The source box can be swivelled to take X-ray Halogen Examination Lamp (12 Volts, 50Watts)
To provide accurate assessment of Infant colour and sufficient illumination during procedures Infant Care Trolley
Transparent acrylic, drop down and lockable side panels Baby Tray with mattress, along with head up/down facility and X-ray cassette holder Mattress density approx 21-25Kg/m³ and with removable, washable, waterproof cover Acrylic X-ray tray with X-ray guide for positioning the cassette 2 Lower shelves Mounted on wheel (100mm) The design aspects also allow easy of dismantling baby care area during cleaning or fumigation
Essential Attachments Side Mounted Monitor tray, Pole for Syringe Pump / Infusion Pump S.S. IV Stand. 4 swivel castors (100mm) - 2 castors with brake



Samco Medical

Website: www.samcomedical.com,