



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**

Infant Radiant Warmer with Detachable Baby Cradle with Microprocessor based temperature controller with 2 modes Skin/Manual mode

**Product Code :**

SAMCOMEDICAL410001A032

**Description :**

Infant Radiant Warmer with Detachable Baby Cradle with Microprocessor based temperature controller with 2 modes Skin/Manual model This is mobile Warmer Stand mounted on wheels with detachable infant care trolley

**Technical Specification**

**Heater Box** Heating Element : Ceramic Heater 650W with parabolic reflector, protected by metallic grid The source box can be swivelled to take X-ray Halogen Examination Lamp (12 Volts, 50Watts)

To provide accurate assessment of Infant colour and sufficient illumination during procedures Infant Care Trolley

Transparent acrylic, drop down and lockable side panels Baby Tray with mattress, along with head up/down facility and X-ray cassette holder Mattress density approx 21-25Kg/m<sup>3</sup> and with removable, washable, waterproof cover Acrylic X-ray tray with X-ray guide for positioning the cassette 2 Lower shelves Mounted on wheel (100mm) The design aspects also allow easy of dismantling baby care area during cleaning or fumigation

**Essential Attachments** Side Mounted Monitor tray, Pole for Syringe Pump / Infusion Pump S.S. IV Stand. 4 swivel castors (100mm) - 2 castors with brake



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),