



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Resuscitation Unit with Fixed Baby Cradle (Standard) with 3 modes Skin / Air / Manual mode

Product Code :

SAMCOMEDICAL410001A041

**Description :**

Resuscitation Unit with Fixed Baby Cradle (Standard) with 3 modes Skin / Air / Manual mode This is mobile Neonatal Resuscitation Unit mounted on wheels with fixed baby cradle & lower shelf along with facility of Oxygenation, Suction and Intubation Mechanical Specification Heater Box Consists of Ceramic Infra Red Heater 650 Watts Metallic Heater Box. The construction of the stainless steel reflector and Guard rail is such that the complete cleaning of heater & reflector is possible without any electrical shock hazards Halogen Examination Lamp (12 Volts, 50 Watts) To provide accurate assessment of Infant colour and sufficient illumination during procedures. Baby Tray (Stainless Steel) Clear acrylic collapsible side panels . Stainless steel baby tray with transparent slot to incorporate the facility of under surface phototherapy. Acrylic tray and cushion mattress with head up/down facility. Acrylic X-ray tray with X-ray guide for positioning the cassette. The design aspects also allow ease of dismantling baby care area during cleaning or fumigation Essential Attachments Lower Shelf S.S. IV Stand. Oxygenation : Oxygen regulator with Hose pipe, Place to keep Cylinder, Chain to hold Oxygen Cylinder, Flow meter, Humidifier bottle, Silicone Patient Tubing, Manual Silicone Resuscitator 250ml with 00 & 01 size mask. Suction : Electrically operated low pressure Suction unit with 2 Autoclavable 500 ml suction jars, Regulator, Suction gauge and silicone tubing. Intubation : Laryngoscope with size 0 & 1 Straight Blades, Silicone Airways 000, 00 & 0



Samco Medical

Website: www.samcomedical.com,

