



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Mouth to Mask Resuscitation

Product Code :
SAMCOMEDICAL410001A043



Description :

Mouth to Mask Resuscitation n ideal solution is to manage PPV in a simplest and safest way. Blow Safe[®] 1/2 is a scientific device designed to deliver PPV with pressure release valve and can be easily used by any Health Care personal. Autoclavable silicone corrugated extension tubing, polycarbonate one-way non-breathing valve, pressure release valve and silicone face mask makes the devices to be used in any ruff-n-tuff condition. Mouth Piece Silicone Corrugated Autoclavable Extension Tubing. The corrugate section prevents blockage of air inlet even when the tubing is in bending position. Due to high volume it creates sufficient pressure than ordinary tubing. Safety Pressure Release Valve One-Way Non-Rebreathing Valve This prevents mixing of inhaled and exhaled air. Silicone Face Mask This provides air tight sealing above nose and mouth. Size: 01



Samco Medical

Website: www.samcomedical.com,