



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Child Resuscitation Kit

Product Code :
SAMCOMEDICAL410001A045



Description :

Child Resuscitation Kit An Essential Carrying Kit for providing Effective Resuscitation to a Child. (Available in a Carry Bag). Contents of the Kit Child Resuscitator, Foot Suction, Laryngoscope, Airways. Child Resuscitator: Rugged 100% Autoclavable; Reusable 500 ml silicone bellow Non rebreathing valve with 40 cm of H₂O pressure release Mask sizes : No. 01 & 02 360° ½ swiveling patient connector with standard 15mm/ 22mm diameter Corrugated PVC oxygen reservoir 1.5 mtr. PVC oxygen tubing Foot Suction: Compact, light weight, easy to handle and operate Durable rubber bellow Long lasting stainless spring to provide friction free pumping Completely autoclavable polycarbonate vacuum jar with lid (500 ml capacity) Scratch proof powder coated frame Laryngoscope : Stainless steel straight blades of size 0 & 1 with penlight handle. Airways : Silicone, autoclavable & reusable in sizes 000, 00 & 0



Samco Medical

Website: www.samcomedical.com,