



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**

Infant Incubator Double Wall (Double Canopy) with drawers with Microprocessor based temperature controller with Skin / Manual mode

**Product Code :**

SAMCOMEDICAL410001A051

**Description :**

Infant Incubator Double Wall (Double Canopy) with drawers with Microprocessor based temperature controller with Skin / Manual mode Technical Specification: Acrylic Doubled Wall Canopy with six port holes. Front Loading. Four Inlets for administration of O2, IV lines or Ventilator tubing. Sliding Baby Tray with External jerk free gear based Head Up / Down facility. X-Ray cassette sliding facility. Stainless Steel rust free Heater, blower sub assembly with changeable micro Air filters. Provision of adequate Storage by 3 drawers and an open shelf. Humidity heater and Hygrometer. Pole for Syringe Pump / Infusion Pump & S.S. IV Stand



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),