



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Infant Incubator Single Wall (Single Canopy) without drawers with Microprocessor based temperature controller with Skin / Air / Manual mode

Product Code :

SAMCOMEDICAL410001A052



Description :

Infant Incubator Single Wall (Single Canopy) without drawers with Microprocessor based temperature controller with Skin / Air / Manual mode Technical Specification: Acrylic Single Wall Canopy with six port holes. Front Loading. Four Inlets for administration of O2, IV lines or Ventilator tubing. Sliding Baby Tray with External jerk free gear based Head Up / Down facility. X-Ray cassette sliding facility. Stainless Steel rust free Heater, blower sub assembly with changeable micro Air filters. Provision of an open shelf in the trolley. Humidity heater and Hygrometer



Samco Medical

Website: www.samcomedical.com,