



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Oxygen Hoods (Medium)

Product Code :
SAMCOMEDICAL410001A054



Description :

Oxygen Hoods (Medium) Single piece, Round Shape, made of autoclavable Polycarbonate material. Trauma free Silicone neck adjustment flap. Bilateral oxygen nozzle (prevents direct flow of cold oxygen on patient's head. Height : 7.5in¹/₂. Width : 8.5



Samco Medical

Website: www.samcomedical.com,