



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Phototherapy Stand

Product Code :
SAMCOMEDICAL410001A055



Description :

Phototherapy Stand Mobile Stand with height adjustment facility along with easy swiveling of source box. Unit can be used along with Infant Care Trolley, Radiant Warmer and Incubator 4 Special Blue CFL Tubes (20W each) 2 Daylight CFL Tubes (20W each) Single Head, surface Size, approx. 0.5m x 0.33m Head Height Adjustment, approx. 0.95 to 1.5m High irradiance at Skin level up to 40uw/cm²/nm (variable with distance) Wavelength - 420 to 500nm, with highest intensity at 470 nm Dual digital cumulative hour timer for tubes usage and patient exposure Antistatic Castors, 2 with brakes Power requirement 1/2 230V/50Hz Power consumption 1/2 100W



Samco Medical

Website: www.samcomedical.com,