



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Phototherapy Stand with Trolley (Single Surface)

Product Code :
SAMCOMEDICAL410001A056



Description :

Phototherapy Stand with Trolley (Single Surface) Technical Specification 4 Special Blue CFL Tubes (20W each) 2 Daylight CFL Tubes (20W each) Single Head, surface Size, approx. 0.5m x 0.33m Head Height Adjustment, approx. 0.95 to 1.5m High irradiance at Skin level up to 40uw/cm²/nm (variable with distance) Wavelength - 420 to 500nm, with highest intensity at 470 nm Dual digital cumulative hour timer for tubes usage and patient exposure Antistatic Castors, 2 with brakes Power requirement $\approx \frac{1}{2}$ 230V/50Hz Power consumption $\approx \frac{1}{2}$ 100W Infant Care Trolley Mobile Bassinet Cot for neonates Cart on 4 swivel castors (100mm Dia), 2 castors with brake Baby Tray with mattress alongwith head up/down facility and X-Ray cassette holder. Mattress density approx. 21-25 Kg/m² and with removable, washable, waterproof cover Transparent acrylic, drop down and lockable side panels Two shelves below baby bed for under surface phototherapy and storage



Samco Medical

Website: www.samcomedical.com,