



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Infant Radiant Warmer (ECO-BASIC)

Product Code :
SAMCOMEDICAL410001A058



Description :

Infant Radiant Warmer (ECO-BASIC) Technical Specification Microprocessor based servo controlled unit allow two operation mode, automatic and manual mode. Audio and Visual alarm function for Power failure, temperature deviation, over temperature & Temperature sensor failure. Transparent Side panel. Heater consists of anti-explosion micro crystal quartz infrared radiation tubes, radiation source 800W. A second thermal cut-out function for more safety. Skin temperature sensor failure function to avoid over temperature



Samco Medical

Website: www.samcomedical.com,