



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Infant Radiant Warmer (ECO-TOP)

Product Code :
SAMCOMEDICAL410001A059



Description :

Infant Radiant Warmer (ECO-TOP) Technical Specification Microprocessor controlled unit allow three operational mode: Pre-warm, Manual or Automatic(Servo Controlled temperature System). Audio and Visual alarm function for Power failure, temperature deviation, over temperature & Temperature sensor failure. Transparent Side panel. Heater consists of anti-explosion micro crystal quartz infrared radiation tubes, radiation source 800W. Heater head can be adjusted $\pm 90^\circ$ horizontally. Convenient X-ray cassette board under infant bed. Inclination of infant bed is adjustable. A second thermal cut-out function for more safety. Skin temperature sensor failure function to avoid over temperature



Samco Medical

Website: www.samcomedical.com,