



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Infant Radiant Warmer (STANDARD)

**Product Code :**  
SAMCOMEDICAL410001A060



**Description :**

Infant Radiant Warmer (STANDARD) Technical Specification Microprocessor controlled unit allow three operational mode: Prewarm, Manual or Automatic(Servo Controlled temperature System). Audio and Visual alarm function for Power failure, temperature deviation, over temperature & Temperature sensor failure. APGAR timer to record therapy time Convenient X-ray cassette under infant bed Heater consists of anti-explosion micro crystal quartz infrared radiation tubes, radiation source 800W. Heater head can be adjusted  $\pm 1/2$  90°  $\pm 1/2$  horizontally Inclination of infant bed is adjustable. Observe light angle is adjustable. Foldable transparent side panels Optional Parts 1. RS232 Connector 2. Drawer 3. Height Adjustable



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),