



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Infant Radiant Warmer (MEDIUM)

Product Code :
SAMCOMEDICAL410001A061



Description :

Infant Radiant Warmer (MEDIUM) Technical Specification : Microprocessor controlled unit allow three operational mode, Pre-warm, Manual or Automatic(Servo Controlled temperature System). Audio and Visual alarm function for Power failure, temperature deviation, over temperature & Temperature sensor failure. APGAR timer to record therapy time Convenient X-ray cassette under infant bed. Heater consists of anti-explosion micro crystal quartz infrared radiation tubes, radiation source 800W. Heater head can be adjusted $\pm 90^\circ$ horizontally. Inclination of infant bed is adjustable. Observe light angle is adjustable. Foldable transparent side panels. LED Phototherapy Unit Optional part: 1. RS232 connector. 2. Drawer



Samco Medical

Website: www.samcomedical.com,