



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Infant Radiant Warmer (PREMIUM)

Product Code :
SAMCOMEDICAL410001A062



Description :

Infant Radiant Warmer (PREMIUM) Technical Specification Microprocessor controlled unit allow three operational mode: Prewarm, Manual or Automatic(Servo Controlled temperature System). Digital skin temperature sensor, the data is exact and stabilizability. Air temperature display, Times count up down. Audio and Visual alarm function for Power failure, temperature deviation, over temperature & Temperature sensor failure. Control panel displayed by VFD screen for precise and stable data. High sensitivity and untouched detector to inspect mattress temperature APGAR timer to record therapy time Observe light angle can be adjusted $\pm 90^\circ$, and the brightness adjusted Heater consists of anti-explosion micro crystal quartz infrared radiation tubes, radiation source 800W. Heater head can be adjusted $\pm 90^\circ$ horizontally Electric infant bed tilt angle is $\pm 10^\circ$ Foldable transparent side panels. Convenient X-ray cassette under infant bed Electrical height adjustment of infant bed. A second thermal cut-out function for more safety Skin temperature sensor failure function to avoid over temperature Optional Part: 1. RS232 connector. 2. Oxygen system. 3. Phototherapy Unit



Samco Medical

Website: www.samcomedical.com,