



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Resuscitation Unit (STANDARD)

Product Code :
SAMCOMEDICAL410001A063



Description :

Resuscitation Unit (STANDARD) Technical Specification: Microprocessor controlled unit allow three operational mode, Pre-warm, Manual or Automatic(Servo Controlled temperature System). Audio and Visual alarm function for Power failure, temperature deviation, over temperature & Temperature sensor failure. APGAR timer to record therapy time Convenient X-ray cassette under infant mattress. Heater consists of anti-explosion micro crystal quartz infrared radiation tubes, radiation source 800W. Heater head can be adjusted $\pm 1/2$ 90° $\pm 1/2$ horizontally Inclination of infant bed is adjustable. Observe light angle is adjustable Foldable transparent side panels. A second thermal cut-out function for more safety. Skin temperature sensor failure protect function avoid over temperature Temperature deviation can be amended in front panel directly Low pressure suction unit, oxygen supply system and LED phototherapy unit Electric height adjusted suit for operation Optional Part: 1. RS232 connector 2. Drawer



Samco Medical

Website: www.samcomedical.com,