



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Resuscitation Unit (PREMIUM)

Product Code :
SAMCOMEDICAL410001A064



Description :

Resuscitation Unit (PREMIUM) Microprocessor controlled unit allow three operational mode, Pre-warm, Manual or Automatic (Servo Controlled temperature System). Audio and Visual alarm function for Power failure, temperature deviation, over temperature & Temperature sensor failure. APGAR timer to record therapy time Convenient X-ray cassette under infant bed Heater consists of anti-explosion micro crystal quartz infrared radiation tubes, radiation source 800W. Heater head can be adjusted $\pm 90^\circ$ horizontally Inclination of infant bed is adjustable. Observe light angle is adjustable Foldable transparent side panels. A second thermal cut-out function for more safety Skin temperature sensor failure function to avoid over temperature. Low pressure suction unit, oxygen supply system and LED phototherapy unit Optional Part: 1. RS232 connector 2. Drawer



Samco Medical

Website: www.samcomedical.com,