



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Infant Incubators with cupboards with baby skin mode and oxygen inlet is optional

Product Code :

SAMCOMEDICAL410001A067

**Description :**

Infant Incubators with cupboards with baby skin mode and oxygen inlet is optional Technical Specification
Microprocessor based servo controlled temperature system Controlled mode: Air Air temperature control range: $25\frac{1}{2}^{\circ}\text{C}$ - $38\frac{1}{2}^{\circ}\text{C}$ ($>37\frac{1}{2}^{\circ}\text{C}$ temperature setting.) Temperature and Heater output displayed by LED Audio and Visual alarm function for Power failure, temperature deviation, over temperature, temperature sensor failure and fan failure. Inclination of infant bed is adjustable. Humidity is continuously adjustable 4 operating windows and 2 access ports. A second thermal cut-off function for more safety Standard Configuration: Main body includes transparent hood, control system, infant bed, water trough, IV pole, Mattress, Support bracket, Castors With optional function Phototherapy up side unit Phototherapy down side unit RS232 Connector Tray & Pole Infant bed pulled from front and side door. Double wall hood



Samco Medical

Website: www.samcomedical.com,