



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Infant Incubators (STANDARD)

Product Code :
SAMCOMEDICAL410001A068



Description :

Infant Incubators (STANDARD) Technical Specification Microprocessor based servo controlled temperature system Controlled mode: Air Air temperature control range : $25\frac{1}{2}^{\circ}\text{C}$ - $38\frac{1}{2}^{\circ}\text{C}$ ($>37\frac{1}{2}^{\circ}\text{C}$ temperature setting.) Temperature and Heater output displayed by LED. Audio and Visual alarm function for Power failure, temperature deviation, over temperature, temperature sensor failure and fan failure. Inclination of infant bed is adjustable. With Phototherapy unit and tray. Humidity is continuously adjustable Double wall hood. 4 operating windows and 2 access ports. A second thermal cut-off function for more safety Standard Configuration: Main body including the transparent hood, control system, infant bed, water trough, IV pole, mattress, cabinets, Silence Castors With optional function Phototherapy up side unit Phototherapy down side unit RS232 Connector. Tray & Pole. Infant bed pulled from front and side door



Samco Medical

Website: www.samcomedical.com,