



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Infant Incubator (LUXURIOUS)

Product Code :
SAMCOMEDICAL410001A069



Description :

Infant Incubator (LUXURIOUS) Technical Specification Microprocessor based servo controlled temperature system Controlled mode: Air and baby skin mode. Air temperature control range: 25°C - 38°C ($>37^{\circ}\text{C}$ temperature setting.) Timer and Heater output displayed by LED Audio and Visual alarm function for Power failure, temperature deviation, over temperature, temperature sensor failure and fan failure. Inclination of infant bed is adjustable. Double wall hood. With phototherapy unit and tray. Humidity is adjustable continuously Standard Configuration Main body (including the transparent hood, control system, infant bed) Skin temperature sensor, Water trough pole, mattress, cabinets with three drawers, Silence Castors With optional function Phototherapy up side unit Phototherapy down side unit RS232 Connector. Tray & Pole. Double wall hood. Infant bed pulled from front and side door



Samco Medical

Website: www.samcomedical.com,