



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Infant Incubator (ADVANCED)

Product Code :
SAMCOMEDICAL410001A071



Description :

Infant Incubator (ADVANCED) Technical Specification Microprocessor based servo controlled temperature system Controlled mode: Air and baby skin mode. Air temperature control range: $25\frac{1}{2}^{\circ}\text{C}$ - $38\frac{1}{2}^{\circ}\text{C}$ ($>37\frac{1}{2}^{\circ}\text{C}$ temperature setting.) Two skin probes for measurement of central and peripheral temperature. Baby weight, X-ray cassette, humidity control system. Temperature, humidity, timer, baby weight and Heat power separately displayed by LED Audio and Visual alarm function for Power failure, temperature deviation, over temperature, temperature sensor failure, baby weight and fan failure. Inclination of infant bed is adjustable. 4 operating windows and 2 access ports. A second thermal cut -off function for more safety With optional function Phototherapy up side unit RS232 Connector



Samco Medical

Website: www.samcomedical.com,