



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Laryngoscope set with Standard Illumination - in Satin/Dull  
Finish

**Product Code :**  
SAMCOMEDICAL7006A008



**Description :**

Laryngoscope set with Standard Illumination - in Satin/Dull Finish The set of four satin blades, includes: Curved macintosh type, size: 1, 2 & 3. Straight miller type, size: 1. Spare bulbs - 4 pcs., & Stainless steel satin handle uses two alkaline batteries (medium size, 1.5v battery)



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),