



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Small Vaccine Carrier, Short Range

Product Code :
SAMCOMEDICAL220001A002



Description :

Small Vaccine Carrier, Short Range Vaccine Storage Capacity as per Std. PQS load : 0.80 Litres. Gross Storage Capacity : 0.80 Litres. Weight Fully Loaded : 2.18 Kgs. Weight Empty (with empty ice pack) : 1.34 Kgs. External surface material : HDPE* Internal lining material : HIPS* Insulation material : CFC free Polyurethane Insulation thickness : 32/36 mm External Dimensions(LxWxH) : 25.0x18.0x21.0 cms. Internal Dimensions(LxWxH) : 17.30x10.30x12.10 cms. Vaccine Storage Dimension : 17.30x10.30x4.50 cms. Lid Type and Fixings : Removable Number of Ice Packs required : 2. Number of Ice Packs supplied : 2. Ice-pack type compliant with PQS E0S/IP01.1 : 0.4 L. Volume of Ice Pack : 0.40 Litres. Min. Cold life Req. : 15 hrs. Conforming to WHO PQS Ref. No. : E04/VC01.1. Cold Life without openings at +43 $\pm$ 1/2 $^{\circ}$ C : 20 hrs.16 mins. Warm Life at -20 $\pm$ 1/2 $^{\circ}$ C : 7 hrs. 55 mins. Cool Life at +43 $\pm$ 1/2 $^{\circ}$ C : 3 hrs. 53 mins



Samco Medical

Website: www.samcomedical.com,