



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Large Vaccine Carrier, Long Range

Product Code :
SAMCOMEDICAL220001A003



Description :

Large Vaccine Carrier, Long Range Vaccine Storage Capacity as per Std. PQS load : 1.35 Litres. Gross Storage Capacity : 1.35 Litres. Weight Fully Loaded : 4.36 Kgs. Weight Empty (with empty ice pack) : 2.22 Kgs. External surface material : HDPE* Internal lining material : HIPS* Insulation material : CFC free Polyurethane Insulation thickness : 35/40 mm External Dimensions(LxWxH) : 24.0x24.0x30.0 cms. Internal Dimensions(LxWxH) : 9.9x9.9x16.50 cms. Vaccine Storage Dimension : 9.0x9.1x16.50 cms. Lid Type and Fixings : Removable Number of Ice Packs required : 4. Number of Ice Packs supplied : 4. Ice-pack type compliant with PQS E0S/IP01.1 : 0.4 L. Volume of Ice Pack : 0.40 Litres. Min. Cold life Req. : 30 hrs. Conforming to WHO PQS Ref. No. : E04/VC01.1. Cold Life without openings at +43 $\frac{1}{2}$ C : 40 hrs.02 mins. Warm Life at -20 $\frac{1}{2}$ C : 15 hrs. 52 mins. Cool Life at +43 $\frac{1}{2}$ C : 7 hrs. 48 mins



Samco Medical

Website: www.samcomedical.com,