



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Large Cold Box Long range

Product Code :
SAMCOMEDICAL220002A002



Description :

Large Cold Box Long range Vaccine Storage Capacity as per Std. PQS load : 18 Litres. Gross Storage Capacity : 22.21 Litres. Weight Fully Loaded : 45.86 Kgs. Weight Empty (with empty ice pack) : 21.78 Kgs. External surface material : LLDPE* Internal lining material : LLDPE* Insulation material : CFC free Polyurethane Insulation thickness : 100 mm External Dimensions(LxWxH) : 70.0x61.0x51.0 cms. Internal Dimensions(LxWxH) : 53.0x37.5x30.5 cms. Vaccine Storage Dimension : 45.5x30.5x16.0 cms. Lid Type and Fixings : Hinged Number of Ice Packs required : 50. Number of Ice Packs supplied : 50. Ice-pack type compliant with PQS E0S/IP01.1 : 0.3 L. Volume of Ice Pack : 0.30 Litres. Min. Cold life Req. : 96 hrs. Conforming to WHO PQS Ref. No. : E04/CB01.1. Cold Life without openings at +43 $\frac{1}{2}$ °C : 126 hrs. 32 mins. Warm Life at -20 $\frac{1}{2}$ °C : 44 hrs. 48 mins. Cool Life at +43 $\frac{1}{2}$ °C : 28 hrs. 05 mins



Samco Medical

Website: www.samcomedical.com,