



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Large Cold Box Short range

Product Code :
SAMCOMEDICAL220002A004



Description :

Large Cold Box Short range Vaccine Storage Capacity as per Std. PQS load : 16 Litres. Gross Storage Capacity : 22.91 Litres. Weight Fully Loaded : 35 Kgs. Weight Empty (with empty ice pack) : 14.50 Kgs. External surface material : LLDPE* Internal lining material : LLDPE* Insulation material : CFC free Polyurethane Insulation thickness : 65 mm External Dimensions(LxWxH) : 64.5x64.0x35.0 cms. Internal Dimensions(LxWxH) : 48.0x47.9x22.6 cms. Vaccine Storage Dimension : 39.9x39.6x14.5 cms. Lid Type and Fixings : Hinged Number of Ice Packs required : 24. Number of Ice Packs supplied : 24. Ice-pack type compliant with PQS E0S/IP01.1 : 0.6 L. Volume of Ice Pack : 0.6 Litres. Min. Cold life Req. : 48 hrs. Conforming to WHO PQS Ref. No. : E04/CB01.1. Cold Life without openings at +43 $\pm$ 1/2 $^{\circ}$ C : 93 hrs. 33 mins. Warm Life at -20 $\pm$ 1/2 $^{\circ}$ C : 31 hrs. 22 mins. Cool Life at +43 $\pm$ 1/2 $^{\circ}$ C : 20 hrs. 42 mins



Samco Medical

Website: www.samcomedical.com,